



CARL ZEISS Oberkochen/Württ.

DEGENHARDT & CO. LTD.


6 CAVENDISH SQUARE, LONDON, W. I.
TELEPHONE NO. LANGHAM 6097/8/9

40-085/1-e

Printed in Germany

Scho. VIII/56 Pooo

ZEISS



INSTRUMENTS FOR DOCTORS AND HOSPITALS

ZEISS - instruments

are aids for the doctor

Ever since Ernst Abbe found his way from the lecturing desk into the work shop of his friend Carl Zeiss, the reciprocal effect of science and practice has stood its test beyond the scope of this enterprise; it has been fused into a symbolic unity under the name of

ZEISS

Scientists and technicians of the Carl-Zeiss works in Oberkochen developed instruments in the fields of microscopy, of optical measuring technics and of medical optics which have become indispensable aids for both doctors and hospitals.

With this publication we present a summary of the most important instruments for doctors and hospitals.



CONTENTS

Microscopes and Accessories

STANDARD JUNIOR Microscope	4
STANDARD Microscope	5
STANDARD MICROMAT	6
ULTRAPHOT II	7
Micro-Projection Equipment	8
Attachment Camera	9
Fluorescence Equipment	10
Stereo-Microscope	11

Optical Measuring Instruments

Percentage Polarimeter	12
Circle Polarimeter 0.05°	12
Extinction Marker II, Electrophoresis Evaluating Instrument	13
Pulfrich Photometer (Pupho), Model G/2	14
Electrophotometer ELKO III	15
Filter Flame Photometer PF 5	16
Spectrophotometer PMQ II	17

Medical-Optical Apparatuses and Instruments

Examination and Operation Microscope	18
Colposcope	19
Slit Lamp Apparatus	20
Ophthalmometer	21
Ophthalmic Diagnostic Set	22
Twin Lamp	23
Spherical Mirror Lamp	23
Reflex-free Ophthalmoscope	24
Magnifiers	24

Index of Special Leaflets	25
---------------------------------	----

We also supply	26
----------------------	----



STANDARD JUNIOR Microscope

The STANDARD JUNIOR serves as a reliable microscope in the routine of medical and biological examinations, courses in microscopy and routine work in laboratories. It is particularly well adapted to the operations required in microscopy and has all the properties of a modern microscope.

Low-set adjustment knobs, co-axial coarse and fine adjustment, **inclined eyepiece tubes** – monocular or binocular, **reliable protection of specimen** – with objectives of higher power through spring mount.

It is possible to use attachments with the microscope for phase contrast and dark-field examinations as well as the micro-projection equipment and the attachment camera. Either an illuminating mirror, a plug-in lamp for connection to the microscope base, or a microscope lamp mounted separate from the microscope serve for illumination.

Leaflet 40-130-e




STANDARD Microscope



The STANDARD Microscope is a large microscope stand for universal application with completely built-in illumination and can be equipped with manifold attachments. It can be adapted for any microscopic work, such as bright and dark field observations, epi-illumination examinations, phase contrast microscopy, photomicrography – from the lowest to the highest magnifications – and for microprojection. An excellently plane field of vision is provided by our new **plane-achromats**.

With the lamp inserted in the microscope base the radiation may be exactly delimited according to Köhler's laws. The inclined binocular tube produces a particularly great brightness without additional magnification.

Leaflet 40-120-e




**STANDARD
MICROMAT**
**Photo-
Microscope**



This apparatus represents a combination of a microscope and a fully automatically operating miniature camera. It serves for microscopic observation and can be used at the same time for photomicrography with black and white or coloured film. The time of exposure is automatically regulated. All customary methods of examination in epi- and transmitted light can be applied.

Leaflet 40-331-e



ULTRAPHOT II
**Camera-
Microscope**




This modern apparatus with its built-in automatic camera is distinguished by its clear design and its easy handling. As to its optical and mechanical performance it stands up against most extraordinary requirements and is adapted for all methods of microscopic inspection. The camera with fully automatic exposure device is especially equipped for plates (9x12 cm).

Leaflet 40-451-e



Micro- Projection Equipment




This equipment was constructed in order to animate scientific lectures, to convey microscopic specimens to a wider circle and to be able to explain individual sections unequivocally.

The equipment is distinguished by great brightness of the projection image and its easy handling. If desired, it can be supplied with automatically regulated feed of carbons for the arc lamp. Every ZEISS-microscope having a free opening in its stand can be attached to the projection base.



Attachment Camera




With the attachment camera any microscope can be turned into an easily operated photomicrographic device. Only few manipulations are required to mount it on the microscope and it is immediately ready to take photographs, since the optical equipment used for observation is also used for taking the photographs. The specimen can be observed while the photographs are being taken. In conjunction with the STANDARD Microscope survey photographs can be taken by using Luminars (slightly magnifying objectives). The camera may be either a miniature camera attachment, one of the usual miniature cameras, such as the CONTAX, EXAKTA, ROBOT, LEICA or a special adapter for plates 6.5x9 cm.

Above all for colour photomicrography we recommend our highly sensitive **exposure meters** with which even with darker specimens excellent results are obtained.



Fluorescence Equipment



The large fluorescence equipment for clinical and diagnostic routine examinations in hematology and parasitology is suitable for **blue-violet and ultra-violet excitation in transmitted and epi-illumination as well as for subjective observation and for photomicrography.**

The equipment is provided with a high efficiency mercury maximum pressure burner OSRAM HBO 200 in

conjunction with suitable light filters (SCHOTT). Microscope and lamp are firmly attached to a base-plate. For photomicrography a miniature camera can be swung over the microscope.

Our **small** equipment is simpler and only intended for BV-excitation. The lamp, placed separately from the microscope, is equipped with an OSRAM HBO 74 mercury maximum pressure burner.



Stereo-Microscope




This very efficient instrument equipped with rapid magnification changer is especially suitable for specimen work and tissue examinations.

The convenient eyepiece of the instrument conveys a non-reversed image with an excellent three-dimensional impression. Its range of magnification is 6-160-fold. Merely by turning a drum another magnification stage is set. The working distance is sufficiently big to allow for all specimen work to be performed. A number of attachments can be supplied, such as lighting fixtures for transmitted and epi-illumination.

Pulfrich Photometer (Pupho)

Model G/2



The **Pulfrich Photometer** (Pupho) Model G/2 represents an improved type of the visual photometer which has been well known for about three decades, especially in medical chemistry. Measuring of the extinction or absorption takes place only in one light path through a restrictor on a magnified, easily read scale. The second restrictor hitherto used was replaced by a grey wedge light absorption device. All cell changing devices (for standard and microcells) can be tempered.

In conjunction with the photo-electric Alternating Light attachment **Wepho** subjective measuring errors and fatigue of the eye are eliminated.


Leaflet 50-501-e

Electrophotometer ELKO III



The standard instrument ELKO II operates with an unusually high degree of measuring accuracy which in various cases is dispensable, if other advantages are emphasized. Hence – at a considerably lower price – the **ELKO III** enjoys the special advantage of accelerated measurement by reading the extinction and absorption values on a ground glass plate with light-mark projection as well as a rapid change from filament to mercury lamp, merely by shifting a lever. Likewise the use of small cells besides standard cells may be of great importance in practice. The well-known ZEISS-S-Filters permit the unrestricted application of the **working instructions "Clinical Photometry"**.


Leaflet 50-633-e



Filter Flame Photometer

If no Zeiss Spectrophotometer is available for use in conjunction with the flame attachment for flame photometric determination, the **Flame Photometer PF 5** successfully fulfils this task. The serum which in the first place is to be examined for Na, K and Ca content is atomized after appropriate dilution, mixed with a fuel gas and lighted. Via a special light filter the light of the flame reaches the receiver, a secondary electron multiplier, and throws a light mark on a ground glass scale via a microprojection system. With the help of a calibration table the scale values shown can be converted into concentration.

The apparatus is fed entirely from the mains; all devices required – with the exception of the high pressure cylinders – are built-in and easily accessible.




Spectrophotometer PMQ II

This instrument represents a further development of the well-known Spectrophotometer M 4 Q. Light-mark reading both in the indicator and in the monochromator is only one of the striking features of the new construction. The reading of the measured value and of the wave length is thus free of errors by parallax.

The recordable spectral region reaches from near infrared (2.5μ) to short-wave ultraviolet ($200 m\mu$). By means of an external lever the receiver is rapidly switched over at choice for the ultraviolet or the visible spectral region.

The instrument can be connected to the 220 V 50 cycles A.C. mains, no batteries or accumulators being necessary. The power consumption for the complete instrument amounts to about 320 VA. With different voltage an intermediate transformer is used.

Further attachments: Chromatogram Attachment (50-657/Chro)
 Flame Attachment (50-675/Flaz)
 Fluorescence Attachment (50-657/Fluor)
 Diffuse Reflexion Device RA 2 (50-657/Re)




Examination and Operation Microscope

The Zeiss Operation Microscope is a binocular microscope on a stand provided with casters. It serves for examining and operating patients who may be either in a sitting or lying position. The working distance amounts to 125 or 200 mm at choice. With the magnification changer five different magnifications from $\times 4$ – $\times 64$ can be set. Illumination and observation paths of rays run approximately parallel to each other. The efficient optical system of the microscope in conjunction with a hinged intensive green filter conveys a particularly plastic and contrasty image. An additional filter renders almost daylight illumination.

A photographic attachment with electronic flash enables photographs to be taken without interrupting observation, also on colour film.

Equipment and accessories according to application.

Leaflet 30-311-e



The Colposcope

The modern operating mechanism of the stand makes it possible to use the colposcope with all types of examination chairs.

The optical and mechanical parts of the Zeiss Colposcope meet the requirements of an up-to-date gynaecological practice.


Working distance: 200 mm.

Range of magnification: $\times 6.3$ to $\times 40$.

Equipment: Colposcope with straight tube, colposcope with inclined tube.

Accessories: Photographic attachment with electronic flash.

Leaflet 30-312-e



Slit Lamp Apparatus

The slit lamp apparatus serves for slit-lamp microscopy of the anterior sections of the eye, of the fundus and of the chamber angle. The magnification changer permits the setting of five different magnifications from x 6 to x 40 without having to exchange eyepieces or objectives.

The photographic attachment allows photographs to be taken of the examination findings.

Equipment: Slit lamp apparatus on instrument table, slit lamp apparatus on mechanical stage, slit lamp apparatus on double table.

Accessories: Hruby lens with fixation device, gonioscopy equipment, anterior chamber depth micrometer, photographic attachment with electronic flash.



Ophthalmometer



With the ophthalmometer the corneal astigmatism and the position of the principal astigmatic meridians are determined. The instrument has great measuring accuracy and a device for directed measurement of the smallest corneal regions. The amount of the corneal astigmatism is automatically calculated and appears on the scale of a reading eyepiece.

Equipment: Ophthalmometer on instrument table, ophthalmometer on mechanical stage, ophthalmometer on double table.

Accessories: Fixation device, contact lens holder.


Ophthalmic Diagnostic Set



The ophthalmic diagnostic set consists of the following instruments:

- Ophthalmoscope,
- Hand Refractometer (Skiascope),
- Slit Lamp Magnifier,
- Diascleral Lamp.

The instruments can be supplied separately or in any desired combination as a set. Since all instruments can be mounted on the same handle, a convenient and rapid exchange during examination is possible.



Leaflet 30-210-e


Twin Lamp

The twin lamp produces a very bright, sharply delimited, adjustable illuminated field. It serves as a lamp for examinations and operations in surgery, gynaecology, ophthalmology, etc.

The twin lamp is mounted on the swivel arm (70 cm or 120 cm long) of the stand, which is provided with casters, and can be swung in any desired direction. Screw-on filters are supplied for special examinations.

Intensity of Illumination:

- 20,000 Lux at 50 cm distance,
- 1,000 Lux at 1.8 m distance.



Spherical Mirror Lamp

Operating lamp for minor surgery and illuminating device for the dressing station.

Equipment:

- Spherical mirror lamp on simple floor stand,
- Spherical mirror lamp with swivel arm 70 cm long,
- Spherical mirror lamp with swivel arm 120 cm long.

Twin lamp with swivel arm
120 cm long

Leaflets see Index p. 25

Reflex-free Ophthalmoscope

after Thorner

For specialist examinations of the fundus, in particular for diagnosing and determining general ailments of the patient.



Magnifiers

Luminous Pocket Magnifier

Practical combination of light source with optical system ideal for magnifier observations. Magnifier heads with 4 different magnification stages at choice.

Illumination of object direct from above.



Folding Magnifiers

Achromatic-aplanatic correction;
Plastic case;

Magnification stages x6, x8, x10 or folding magnifier x3 and x6 = x9.



Binocular Head Magnifier

Stereoscopic viewing with x2.25 magnification;
big field of view;
both hands free for working.

Index of Special Leaflets

(listed alphabetically)

In this master leaflet only the most important characteristics of the apparatuses and instruments are enumerated. For detailed information we shall be glad to supply the following leaflets:

Apparatus/Instrument	Leaflet No.
Alternating Light Attachment (Wepho)	50-651-e
Attachment Camera	40-410-e
Binocular Head Magnifier	ZW-108
Circular Polarimeter 0.05°	ZW-360
Colposcope	30-312-e
Electrophoresis Evaluating Instrument	50-610-e
Electrophotometer Elko III	50-633-e
Electrophotometer Elko II	50-631-e
Examination and Operation Microscope	30-311-e
Extinction Marker II	50-610-e
Filter Flame Photometer PF 5	50-620-e
Fluorescence Equipment	40-215-e
Folding Magnifier	ZW-110
Luminous Pocket Magnifier	ZW-112
Micro-Projection Equipment	40-370-e
Ophthalmic Diagnostic Set	30-210-e
Ophthalmometer	30-230-e
Percentage Polarimeter	ZW-350
Pulfrich Photometer (Pupho)	50-501-e
Reflex-free Ophthalmoscope	ZW-102
Slit Lamp Apparatus	30-245-e
Spectrophotometer PMQ II	50-657-e
Spherical Mirror Lamp	30-520-e
STANDARD JUNIOR Microscope	40-130-e
STANDARD Micromat	40-331-e
STANDARD Microscope	40-120-e
Stereo-Microscope	40-710-e
Twin Lamp	30-515-e
Ultraphot II	40-451-e



Leaflets see Index p. 25

We also supply:



Microscopes:

- Travelling Microscopes
- Polarizing Microscopes
- Metal Microscopes
- Plankton Microscopes
- Trichinoscopes

Microscope Accessories:

- Microscope Lamps
- Drawing Apparatuses
- Photo-micrographic Equipments
- Micro-projection Apparatuses
- Micro-cine-Cameras
- Grinding and Polishing Machines
- Heating Stages
- Polarizing Equipments

Ophthalmic and Medical Instruments:

- Eye Microscopes
- Gonioscopes
- Eye Cameras
- Ophthalmoscopes
- Skiascopes
- Slit Lamp Magnifiers
- Diascleral Lamps
- Vertex Refractionometer
- Hand Lamps
- Otoscopes
- Dermatoscopes
- Technoscopes
- Capillary Microscopes
- Polarizing Spectacles

Optical Measuring Instruments:

- Refractometers
- Spectroscopes
- Electrophotometers for Transmission, Reflexion and Turbidity Measurements
- Monochromators
- Registering Spectrophotometer
- Galvanometer Reading Device
- Photoelectric Diffuse Reflexion Photometers
- Hydrogen Lamps
- Spark Generator and Interruption-arc Instrument

Precision Measuring Instruments:

- Universal Measuring Microscopes
- Big Toolmaker's Microscopes
- Small Toolmaker's Microscopes
- Measuring Microscopes 0-50
- Light Section Microscopes
- Big Optical Dividing Heads
- Optical Rotary Tables
- Microscopes for Use on Machine Tools
- Lead Measuring Instruments
- Optical Flat Gauges

Interferometric Instruments:

- Firedamp Interferometers
- Gauge Block Interferometers
- Diffusion Interferometers
- Interference Microscopes
- Interference Surface Tester
- Schlieren Instruments

Surveying Instruments:


- Levels and Theodolites

Spectacle Lenses – Polarizing High Efficiency Objectives – Magnifiers – Polarizing Filters – Colour Filters – Binoculars – Astronomical Instruments – Planetaria – Electron Microscopes

The illustrations are not binding in all details for
the construction of the instruments.

Electros, or reductions of same when available,
are supplied for scientific publications.

Reproduction of illustrations or text without our
consent is not permitted.



CARL ZEISS · Oberkochen/Württ.