hund
WETZLAR

V 100

V 200

V 300

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Instruction for use

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Manual de instrucciones

Istruzioni d'uso

Deutsch

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English

Congratulations for deciding on a quality product from **hund**.

Please compare the contents of this consignment with the enclosed delivery note.

To ensure our factory's guarantee (lifetime for mechanical and optical parts, 12 months for electrical and electronic parts), please send the enclosed guarantee card duly filled in and signed as soon as possible to our address.

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HELMUT HUND GMBH
D - 35580 WETZLAR

has established and applies
a quality system for
Microscopy, Measuring Technologies,
Electronics and Plastics Technology,
Mechanics and Optics

An audit was performed, Report No. 94/1249
Proof has been furnished that the requirements according to
DIN EN ISO 9001
are fulfilled.

The certificate is valid until
June 1998
Certificate Registration No.
73 100 132

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Components of the microscope

Illustration 1 shows the Vb 365 Wilocyt, a microscope in the V series. The items illustrated show the parts that can differ slightly according to the model.

- 1.1 Eyepieces
- 1.2 Binocular tube
- 1.3 Stand
- 1.4 Objectives
- 1.5 Specimen holder
- 1.6 Mechanical stage
- 1.7 Knobs for movement of the mechanical stage in x- and y-direction
- 1.8 Condenser
- 1.9 Mains switch with variable control
- 1.10 Base
- 1.11 Field diaphragm
- 1.12 Knob for coarse and fine movement of the mechanical stage in z-direction
- 1.13 Knob for movement of the condenser
- 1.14 Lever for aperture diaphragm
- 1.15 Revolving nosepiece

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Illustration 2

- 2.1 Eyepiece
- 2.2 Scale for interpupillary distance
- 2.3 Eyepiece tube
- 2.4 Stand
- 2.5 Revolving nosepiece
- 2.6 Objective
- 2.7 Knurled screw
- 2.8 Gripping ledge

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Assembling the microscope

Unpack the stand and place it on a stable, vibration free surface, away from strong light sources. Attach the binocular tube to the stand (ill. 2) securing it with the tube fastening screw (2.7). Remove the protective covers from the eyepiece tubes (2.3) and insert the eyepieces (2.1) into the eyepiece tubes.

Objectives (2.6), mechanical stage (1.6) and condenser (1.9) have already been fixed to the stand (1.3) in our factory. Take all accessories out of the packing and store these in a safe place. Connect the instrument to a plug socket with the mains cable. The microscope is now ready for use.

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Adjusting the microscopic image

– Quick reference guide –

(For detailed information please see the following pages).

- a) Switch on the instrument.
- b) Adjust the binocular tube (ill. 2) (interpupillary distance, eyepiece tubes)
- c) Move the condenser to its upper position, using the appropriate knob (1.13).
- d) Open the diaphragm and the aperture diaphragm (1.11 and 1.14).
- e) Move the mechanical stage (1.6) downwards by means for the coarse movement for the coarse movement for the mechanical stage (1.12).
- f) Place a specimen on the mechanical stage and fix it with the specimen holder.
- g) Rotate the nosepiece (1.15) so that the objective (1.4) with the lowest magnification is in use.
- h) Adjust the light intensity.

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Setting up the microscope

– Detailed instructions – (For quick reference guide see page 8 and 9)

a)

For the V-series microscopes the following methods of illumination are available.

V 100-series

Plano-concave mirror or attachable light source 230 V/15 W for external illumination.

V 200-series

Incorporated illumination 220 V/15 W with ON/OFF switch.

V 300-series

Halogen illumination 6 V/10 W with knob for ON/OFF switch and continuous regulation of the brightness control.

b)

All kinds of tubes, i.e. monocular tube, binocular tube or binocular tube with vertical phototube are fixed to the microscope stand by means of the tube fastening screw. The tube is 360° rotatable and can be adjusted by loosening the tube fastening screw. The binocular tubes have to be adjusted according to the interpupillary distance of the observer. Correct the tube length by transferring the index of the interpupillary distance scale (2.2) to the scale of the eyepiece tubes (2.3). (For the MEDICUS, only the left eyepiece tube can be adjusted).

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c)

The condenser (1.8) can be moved up or down either by means of the knob for the adjustment of the condenser (1.3), or by turning the condenser in the sliding sleeve.

If a Wilocyt condenser (ill.6) is fixed to the instrument turn the revolving disc (6.3) until the white spot is visible and clicks into place under the white marking. This is the position for brightfield.

d)

Open the aperture diaphragm (3.2, 4.2, 5.2, 6.4) by turning the lever on the left. Open the field diaphragm (only for models V 300, see ill., 1.11) by turning the black knurled wheel.

k)

Adjust the interpupillary distance of the eyepiece tubes (2.3) as in b). Focus the microscopic image for the right eye by means of coarse and fine movement (1.12) and adjust the ultimate sharpness for the left eye by turning the eyepiece tube.

l)

Adjusting the illumination method:

- 1) brightfield illumination (all models)
- 2) darkfield illumination (Wilocyt-condenser)
- 3) Phasecontrast (Wilocyt, Wilomed, MEDICUS Ph)

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1) Brightfield (for all models)

Adjust the V 300 with the 10x objective. Close the field diaphragm (1.11). If the condenser mounted is as illustration 5, swing out the lens (5.4) to the left side. Lower the condenser from the upper position until the rim of the field diaphragm is clearly visible. Open the field diaphragm slightly past the field of view. Then close the aperture diaphragm (5.2, 3.2, 4.2, 6.4), until optimal contrast of the image is obtained (approx. 2/3).

2) Darkfield

With the Wilocyt condenser (ill. 6) darkfield illumination can be achieved up to a maximum objective magnification of 40x. With this objective a darkfield drop-in diaphragm has to be inserted. To fit diaphragm, detach the objective from the revolving nosepiece by turning the objective to the left. Rotate the revolving disc (6.3) until the violet spot clicks into position under the white marking (6.2).

Aperture diaphragm and field diaphragm need to be completely opened.

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3) Phasecontrast

For phase contrast illumination two components are always necessary. Special phase contrast objectives (marked Ph) and special condensers (with a phase ring).

Aperture and field diaphragm have to be completely opened.

Illustration 3

For the Wilomed and MEDICUS Ph models, phase contrast is only possible with the objective A Ph 40x. The condenser diaphragm for phase contrast is inserted in the condenser from below. Open the aperture diaphragm (3.2) first.

Illustration 6

For the Wilocyt, phase contrast can be achieved with objectives A Ph 20x and A Ph 40x. The revolving disc should be rotated until the green-blue spot clicks into place under the white marking.

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Hints for maintenance

To ensure your instruments' long service it should be checked and serviced regularly. In case of malfunction, please contact your local dealer or the manufacturer directly.

Some useful hints for your microscope:

No light

- Check the fuse in the plug.
- Check the mains cable and the electrical circuit.

Poor quality image:

- Adjust focus
- Adjust condenser
- Reset instrument as in quick reference guide
- Clean the eyepieces and the objectives (front lenses) with cotton bud and alcohol. Always clean oil objectives thoroughly, immediately after use.

Changing the lamp

- Pull out mains cable
- Remove cover plate underneath the base plate.
- Pull out the lamp (V 300).
- Unscrew the lamp (V 200).
- Insert a new lamp.
- Do not touch the glass body with your fingers

Attention: Lamp failure
Before changing the lamp, wait approx. 5 minutes until the lamp has cooled down, otherwise there is a danger of getting burned.

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Adjusting the instrument V 300

- Place the specimen on the mechanical stage and focus the specimen
- Switch to the 10x objective
- Close the field diaphragm (1.10) to the side.
- Centre the field diaphragm mechanically by moving the knurled wheel (1.10) to the side.
- Move the condenser downwards as in 3c, until the diaphragm is sharply visible.
- Centre with centering keys (Ill.4).
- The model V 200 has been adjusted in the factory. No further adjustment is required.

Practical hints

- Use only the original type of microscope lamps (power and shape).
- Put the specimen on the mechanical stage with the cover glass on top.
- After use cover the microscope with the dust cover.
- Choose a suitable place for your instrument.
- Insert the filters (blue filter, green filter or ground glass filter) in the filterholder on the condenser or in the opening for the light exit.
- Do not use more than one filter at the same time.
- Only use immersion oil for the corresponding objectives (marked „oil“ on the objectives).
- Check and adjust the instrument only when the specimen is in focus.

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Accessories

We offer a wide range of accessories for various fields of application, to fit our instruments.

Should you have any question please do not hesitate to contact us, then we would be pleased to send you information on any of the following:

- image processing
- image analysis
- incident light microscopy
- fluorescence microscopy
- polarisation
- darkfield microscopy
- phase contrast microscopy
- stereomicroscopy
- inverted microscopes
- measuring with the microscope
- microscopes with video cameras + monitor
- documentation
- with SLR-cameras
- with Polaroid cameras
- with PAC
- with video printer
- Accessories
- Eyepieces
- Objectives
- Heating stages
- Attachable mechanical stages
- Micromanipulators
- Cold light sources
- Light guides

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Technical data

Power supply:
230 V – 50-60 Hz
115 V – 50-60 Hz
Power consumption:
max. 15 W
Lamp for V 200: 230 V/15 W
115 V/15 W
Lamp for V 300:
6 V/10 W

Fuses:
(5 x 20 mm n. IEC 127-2/5):
T 0,2 A for 230 V, T 0,4 A for
115 V
Safety standard:
IEC 1010-1/EN61010-1/
VDE 0411 safety class 1
Weight: 5,5 kg
Dimensions:
Height: 38 cm,
Breadth: 17 cm,
Depth: 24 cm

Ambient conditions:

Range of temperature:
5°C.–40°C.
Humidity:
max. 80% for 31°C, Linear
decreasing to 50% for
40°C.
Mains fluctuations:
max +/- 10%.
Contamination level:
2 according to IEC 664.
Transient overvoltage:
Category II.

Safety information:
The microscope corres-
ponds to the safety stan-
dards for electrical meas-
uring and laboratory instru-
ments according to IEC
1010-1/EN 61010-1 and left
the company in technically
sound condition.

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To ensure the proper function and safety of the instrument please note the guidelines in the instruction manual. The microscope must only be used as directed in the instruction manual.

Maintenance and repair work should only be carried out by authorised personnel or our service department. Only original spare parts and components must be used as replacement parts. In case of doubt please contact our representative. No responsibility is taken for failure to comply with the safety guidelines.

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Helmut Hund GmbH